

- Kills SARS-CoV-2 on hard non-porous surfaces
- Kills the virus that causes common cold^{††}
- Kills H1N1 Influenza A Virus
- Ready-to-use disinfectant spray
- Kills 99.9% of germs*

DESCRIPTION:

This product is a one-step*,¹ hospital-use germicidal disinfectant cleaner and deodorant designed for general cleaning, disinfecting, and controlling non-pathogenic mold and mildew on hard, non-porous inanimate surfaces. Quickly removes dirt, grime, non-pathogenic mold and mildew, food residue, blood, and other organic matter commonly found in hospitals. It also eliminates odors leaving non-porous surfaces smelling clean and fresh. Use where odors are a problem.

*When used according to disinfection directions.

¹ “One Step” claim is limited to association with non-food contact hard non-porous surface use.

COOLING COILS DISINFECTANT CLEANER: This product is a ready-to-use, non-abrasive disinfectant cleaner for small air conditioner, commercial air handling unit and HVAC cooling coils and refrigeration equipment evaporation coils. When using this product to disinfect and clean cooling coils and evaporation coils, disconnect power to the unit and allow hard non-porous surfaces to come to room temperature before applying product.

Follow directions under **DISINFECTION, DEODORIZING AND CLEANING** section.

For plastic and painted surfaces, spot test on an inconspicuous area before use.

Evap-Fresh™ No Rinse Evaporator Coil Cleaner and Disinfectant eliminates 99.9%* of SARS-CoV-2 the Coronavirus that causes COVID-19 on hard non-porous surfaces in 60 seconds.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. This product is not for use on critical or semi-critical medical device surfaces. Fill or refill the 32oz. Evap-Fresh™ No Rinse Evaporator Coil Cleaner & Disinfectant spray bottle with Evap-Fresh™ No Rinse Evaporator Coil Cleaner & Disinfectant refill. Pour product into spray bottle over sink or tub basin.

Indoor Air Quality Products

Evap-Fresh™ No Rinse Evaporator Coil Cleaner and Disinfectant



Packaging

| | |
|-----------|---------|
| 1 gallon | 4166-08 |
| 1 quart | 4166-24 |
| 55 gallon | 4166-01 |

DISINFECTION, DEODORIZING AND CLEANING:

Remove visible soil prior to application of the product. Hold container six to eight inches from surface to be treated. Spray area until it is covered with solution. Allow the surface to remain wet for **3 minutes**. No scrubbing is necessary. Wipe off with a clean cloth or allow to air dry. The product will not leave grit or soap scum.

Bactericidal Activity: When used as directed, this product kills the following organisms on hard, non-porous inanimate surfaces when the treated surface is allowed to remain visibly wet for **3 minutes**.

BACTERICIDAL ACTIVITY (BACTERIA): Using AOAC test methods in the presence of 5% blood serum and a 3 minute contact time, kills the following organisms on hard non-porous inanimate surfaces:

| Organism | ATCC Number | Contact Time |
|---|--------------|--------------|
| <i>Pseudomonas aeruginosa</i> | (ATCC 15442) | 3 minutes |
| <i>Staphylococcus aureus</i> | (ATCC 6538) | 3 minutes |
| <i>Salmonella (choleraesuis) enterica</i> | (ATCC 10708) | 3 minutes |
| <i>Escherichia coli</i> | (ATCC 11229) | 3 minutes |
| <i>Streptococcus pyogenes</i> | | 3 minutes |
| <i>Listeria monocytogenes</i> | (ATCC 35152) | 3 minutes |
| <i>Yersinia enterocolitica</i> | (ATCC 23715) | 3 minutes |
| <i>Salmonella (typhi) enterica</i> | (ATCC 6539) | 3 minutes |

Tuberculocidal Activity: This product exhibits disinfectant efficacy against Mycobacterium tuberculosis BCG at 20°C when the treated surface is allowed to remain wet for **3 minutes**.

†Virucidal Activity: When used as directed, this product kills the following organisms on hard, non-porous inanimate surfaces: Human Coronavirus (ATCC VR-740), SARS associated Coronavirus, Avian Influenza A Virus strain H3N2 (ATCC VR-2072), and Avian Influenza A Virus strain H9N2, when the treated surface is allowed to remain visibly wet for **2 minutes**. Rhinovirus type 39 (ATCC VR-340) and Rotavirus (ATCC VR-2018), when the treated surface is allowed to remain visibly wet for **3 minutes**. Hepatitis A Virus (HAV), when the treated surface is allowed to remain visibly wet for **10 minutes**.

Mildewstat: To control non-pathogenic mold and mildew on pre-cleaned, hard, non-porous surfaces spray surface to be treated making sure to wet completely. Let air dry. Repeat application at weekly intervals or when mildew growth appears.

The SDS is available on the Nu-Calgon website at www.nucalgon.com.